AS - IOL

Lente intraocular asférica AS-IOL

AS-IOL aspheric intraocular lens

Asphärische Intraokularlinse

AS-IOL

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AS-IOL ASPHERIC INTRAOCULAR LENS

Description

AS-IOL aspheric intraocular lenses are implants that are indicated for lens replacement in the surgical correction of aphakia.

AS-IOL aspheric lenses are designed to correct the spherical aberration of the cornea.

It is a design which provides optimised optical quality by compensating for the positive spherical aberration of the cornea.

Characteristics

Monoblock aspheric lens.

Reference: LLASF60 - LLASY60.

Optical/haptic material:

LLASF60:

- Hydrophilic aspheric intraocular lens with UV filter and natural filter for blue light (370nm of transmission UV at 10% cut off).
- Haptics in the shape of a "C".
- Angulation 0°.
- Optical diameter: 6.0 mm

- Total length: 12.50 mm
- Strength: -5 +40 (increase +1 dioptre).
+12 +25 (increase +0.5 dioptre).

LLASY60:

- Hydrophilic aspheric intraocular lens with UV filter and natural filter for blue light (370nm of transmission UV at 10% cut off).
- Haptics in the form of a perforated plate.
- Angulation: 0°
- Optical diameter: 6.0 mm
- Total length: 10.75 mm
- Strength: -5 +40 (increase +1 dioptre).
+12 +25 (increase +0.5 dioptre).

The lens is supplied, placed in a holder, inside a glass vial which is sealed with a plug for easy opening and is protected by sterile, individual blister packaging.

The AS-IOL intraocular lenses are sterilised by steam.

Precautions for use

To prevent the IOL from becoming temporarily opaque before being implanted or used, we advise allowing it to achieve a stable temperature of

25°C or higher, as with viscoelastic material, saline solution and medicines.

During anterior segment surgery the following precautions for use are recommended:

- Check the integrity of the individual sterilization protection of the product before using. Do not use if the primary packaging is opened or damaged.
- The AS-IOL intraocular lenses are for single use and should not be re-sterilised under any circumstances by any method.

Known risks of a possible reuse of the product are:

- Reusing this product creates a potential risk of infection for the patient. Contamination of the product may cause injury to the patient.
- Reprocessing of the product may compromise its structural integrity, as well as its sterility and functionality.

The concomitant medications requiring precautions for use are:

- Local and/or general carbonic anhydrase inhibitors because of the antisecretory effect on the ciliary body which reduces the secretion of aqueous humour with a loss of secretion of 50%.
- Broad-spectrum antibiotics (VANCOMYCIN) used in intraocular irrigating solutions which modify the pH and the osmolarity of the aqueous humour.

In case of difficulty when opening the vial, do not force the vial open. Forcing a vial open may break the vial creating sharp edges. Report the incident to the manufacturer..

The eye surgeons are solely responsible for their own surgical techniques.

Incompatibility

Do not use quaternary ammonium (benzalkonium chloride) with AS-IOL intraocular lenses.

Contraindications

Apart from non-specific contraindications related to any form of ocular surgery, the following non-exhaustive list of specific contraindications must be respected:

- Bilateral congenital cataracts.
- Recurrent inflammations of the anterior chamber and/or the posterior chamber.
- Pre-existing chronic glaucoma.
- Microphthalmia.
- Corneal dystrophies, particularly endothelial.
- Diabetic or hypertensive retinopathy.

In principle, patients with negative or no spherical aberration would not benefit (optimisation of optical quality) from the isoplanatic spherical design, since it would increase the negative spherical aberration. Patients with high astigmatism, or with corneas with highly asymmetric elevation maps, would not benefit from the design.

However, there are no contraindications for the use of AS -IOL intraocular lenses provided they are used according to the product information brochure.

Complications related to surgery

Any pre-existing eye conditions may be aggravated by the implantation of an intraocular lens.

Intraoperative complications

- Glaucoma
- Vitreous herniation
- Secondary membrane
- Retrolental membrane
- Retinal detachment
- Iridial atrophy
- Severe ametropia and aniseikonia
- Replacement or removal of the lens
- Excessive intraoperative vitreous loss
- Cystoid macular edema

Postoperative complications

- Corneal edema
- Corneal dystrophy
- Pupillary block
- Iridocyclitis and hyalitis
- Precipitations on the lens
- Cases of corneal opacification of the lens due to calcium precipitates have been reported in patients with preexisting conditions or with managed diseases, such as Diabetes.
- Off-centre lens
- Dislocation and subluxation of the lens
- Endophthalmitis and panophthalmitis
- Hemorrhage

Instructions for use

Implantation technique for AS-IOL intraocular lenses:

- Check the expiry date and the strength of the lens in dioptres.
- Pick up the lens and fold it with the forceps, making sure that the fold is perfectly symmetrical.
 - Rotate with the forceps and irrigate the lens with a sterile solution.
 - Introduce the lens horizontally through the incision.

Model LLASY60: Put the lens in place taking care that the notched haptics maintain the “Z” shape.

Model LLASF60: Put the lens in place taking care that the haptics maintain the “Z” shape.

- Straighten the forceps in the vertical plane.
- Release the lens slowly and withdraw the forceps.

- The lens should be fully unfolded and centred in the capsular bag.
- Use of an injector system. This lens can be implanted using an injector system (e.g. MEDICEL injector or similar), always following the instructions for use of such systems. In this case the maximum time folded may not exceed three minutes.

For professional use only.

Reporting customer complaints including quality complaints, adverse events and other medical device related observations:

Customer complaints including quality complaints, adverse events and other medical device related observations should be reported to AJL without delay. A report describing the details of the complaint/event, the applied therapy, the product type, LOT/serial number of the medical device used is requested.

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AJL OPHTHALMIC, S.A.